



# PDC Color Coded Alert Bands® & Ident-Alert™ Snaps



Hospital standardization is spreading for the use of Color Coded Indicators to warn healthcare staff about patients at risk. Precision Dynamics Corporation's (PDC) **Alert Band® and Ident-Alert™ Colored Snaps** provide easy, accurate, and cost effective patient classification, while meeting state standardization guidelines.

## PDC Alert Bands®

- Extra large text provides classification to clinicians (assisting in dim light conditions or for those who are color blind)
- Only primary and secondary colors used as indicators
- Patented PDC SecurSnap® closure prevents transfer of wristbands
- Custom logos and cues communicate less common conditions
- Soft vinyl & plastic materials are latex-free for optimum patient comfort and safety
- Custom imprinting available upon request



**Alert Band® and Ident-Alert™ Colored Snaps are manufactured in North America and are ISO-9001 certified, reflecting PDC's commitment to ensure all products are designed, constructed, and tested to the highest quality standards.**

## Alert Bands® Ordering Guide

### ClearImage™ Alert Bands®

Info Area: 3" L x 1/2" H 500 Bands/Box


Constructed of a unique translucent material for a distinctive look and feel



### Sentry® Superband® Alert Bands® Narrow

Info Area: 3" L x 1/2" H 500 Bands/Box

Constructed of lightweight non-stretch polyester





### Other Risk Reduction Strategies

- Remove wristbands that have been applied at another facility.
- Initiate banding on admission, changes in condition, and when pertinent information is obtained during the Hospital stay.
- Educate patients and family members about the purpose of color coding.
- Coordinate chart/whiteboard, care plan, door signage, etc. with the same color coding.
- Verify color-coded wristbands on assessment during hand-off of care and during facility-to-facility transfer.
- Use the PDC Alert Band® or Ident-Alert™ Snaps as an "alert". Use the patient, medical record and other resources to verify the intent of the wristband.
- Hospitals recommend that people wearing "social cause wristbands," such as yellow "Live Strong" bracelets or pink breast cancer awareness bracelets remove them while in the hospital to avoid possible confusion among patient wristbands used for healthcare purposes.

## Sentry® Superband® Alert Bands® Adult/Pediatric

Info Area: 2 5/8" L x 1" H 500 Bands/Box

Constructed of lightweight non-stretch polyester



## Sentry® Superband® Alert Bands® Adult

Info Area: 3 7/8" L x 1" H 500 Bands/Box

Constructed of lightweight non-stretch polyester



## Ident-Alert™ Colored Snaps

- Conveniently classify multiple patient conditions—all on one wristband
- Smooth rounded edges for added patient comfort
- Use in conjunction with all PDC snap closure wristbands
- Six standard colors; 100 snaps / bag



### Ident-Alert™ Colored Snaps Ordering Guide

<b>Part #</b>	8615-12-PDG	8615-13-PDG	8615-14-PDG	8615-16-PDG	8615-22-PDG	8615-33-PDG
<b>Color</b>	Pink	Blue	Yellow	Red	Green	Purple
<b>Snaps/Bag</b>	100	100	100	100	100	100



### Precision Dynamics Corporation

Tel 800.772.1122 | 818.897.1111

Fax 818.686.9314

13880 Del Sur Street

San Fernando, CA 91340

Email [info@pdcorp.com](mailto:info@pdcorp.com)

[www.pdcorp.com/healthcare](http://www.pdcorp.com/healthcare)

See possible color coding initiatives in your state and learn more about patient safety through color alerts.

Visit [www.alertbands.com](http://www.alertbands.com)

Call **800.772.1122**