



Identification Solutions for Blood Banks



POSITIVE ID. POSITIVE OUTCOMES.™



ID SOLUTIONS FOR BLOOD BANKS

With more than 40,000 units of blood needed each and every day, blood banks play a critical role in saving lives and helping hospitals deliver quality patient care. There is no room for error or time to slow down. With that in mind, Precision Dynamics offers products that deliver positive identification and increased efficiency. We understand the importance of positive identification in blood banks and have been providing Positive ID and Positive Outcomes™ for over 50 years.




SECURLINE® BAR CODE BLOOD BAND

Meets The Joint Commission's National Patient Safety Goal for blood transfusions for a one-person verification system


- Streamlines the blood transfusion process and improves patient safety by reducing human errors.
- Saves time and money by streamlining the patient ID and blood transfusion workflow and minimizing paperwork.
- Works with your existing blood bank system, so there's no IT involvement or capital investment.
- Helps prevent life threatening and costly misidentification errors.
- Choose from Write-On or Label Styles.

LABEL STYLE 2410-**-PDB



- Designed for **inpatient or outpatient use**
- Accommodates pre-printed labels up to 3" L x 1" H
- Clear adhesive shield on the wristband's information area protects patient ID information and bar code from solvents for accurate bar code scanning
- Total of 12 matching labels (128 bar code symbology, alpha-numeric serialization) (8 on draw tube tail and 4 on wristband)
- Advanced material resists moisture, alcohol and hand sanitizer to preserve bar code integrity for accurate scanning

WRITE ON STYLE 2420-**-PDB



- Designed for **inpatient use**
- Allows you to write patient information one time only on the draw tube label while the non carbon paper transfers patient information onto the wristband
- Total of 14 matching labels (128 bar code symbology with alpha-numeric serialization) (8 on draw tube tail, 5 on wristband, and 1 on draw tube label)

PRODUCT NUMBER	SIZE OF BAND	INFORMATION AREA	MAX LABEL SIZE	BANDS PER BOX
2410-**-PDB	12" L	3 1/4" L x 1 1/4" H	3" L x 1" H	150 bands (15 dispenser boxes of 10 bands each)
2420-**-PDB	12" L	2" L x 1 1/4" H	N/A	150 bands (15 dispenser boxes of 10 bands each)

Substitute ** for numbered color: 14 Yellow 22 Kelly Green 16 Red 17 Orange 18 Light Blue 26 Gray

**POSITIVE ID.
POSITIVE OUTCOMES.**

CUSTOMER SUCCESS STORY

"The decision to switch to Securline Bar Code Blood Band was easy. We wanted a band that was easy to apply, secure, and comfortable and that would work with our existing bar code transfusion system. The Securline product made our jobs easier while also giving our patients peace of mind and comfort."

-Mary Ann Sharpe
Blood Bank Manager
Research Medical Center,
HCA Midwest Health System
Kansas City, MO

SECURLINE® BLOOD BANDS

The top wristband choice of hospitals and blood centers to reduce misidentification errors, while increasing patient comfort and safety.

- Patented SecurSnap® closure prevents transfer and tampering
- Easy-to-Read Data Area: a large, white surface offers clear and easy readability which can help reduce medical errors
- Efficient Easy-Peel Labels: includes 23 easy-peel labels (more than any other blood band style) for quick ID and tracking of specimens
- Each box includes one Xtender Band™ for larger patients
- Intended for in-patient use
- For outpatient use, please use Blood Band with Safeguard® Sealident® Super Wide 881

SECURLINE® BLOOD BAND



PRODUCT NUMBER	SIZE OF BAND	DESCRIPTION	BANDS PER BOX
2403-**-PDB	17" L	Securline Blood Band	150 bands (15 dispenser boxes of 10 bands each)

Substitute ** for numbered color:



SECURLINE® BLOOD BAND PLUS

- Features adhesive peel-off tail that applies directly to draw tubes



PRODUCT NUMBER	SIZE OF BAND	DESCRIPTION	BANDS PER BOX
2406-**-PDB	17" L	Securline Blood Band Plus	150 bands (15 dispenser boxes of 10 bands each)

- New DOB field is an additional identifier

new!



PRODUCT NUMBER	SIZE OF BAND	DESCRIPTION	BANDS PER BOX
2407-**-PDB	17" L	Securline Blood Band Plus	150 bands (15 dispenser boxes of 10 bands each)

Substitute ** for numbered color:



**POSITIVE ID.
POSITIVE OUTCOMES.**

**CUSTOMER SUCCESS
STORY**

“Our work is based on having the correct blood sample for the right patient. We can’t afford not to use a reliable blood band as the integrity of our work and the safety of our patients depend on it.”

- Linda Harvey
Blood Bank Supervisor

Northeast Alabama Regional
Medical Center
Anniston, Alabama




PDC'S BLOOD BANDS MEET THESE JOINT COMMISSION PATIENT SAFETY GOALS & STANDARDS:

- Goal #1: Improve the accuracy of patient identification
- Goal #2: Improve the effectiveness of communication among caregivers
- Goals #13: Encourage patients' active involvement in their own care as a patient safety strategy

ALL PRODUCTS ARE LATEX-FREE

SECURLINK® BLOOD RECIPIENT SYSTEM

The Securlink® System links the patient, the transfused blood and all the related paperwork to assure safe and complete blood transfusion. Order the full system (includes 250 bands and forms) or separately. The wristband is made of soft vinyl, featuring a tamper proof patented SecurSnap® closure.



PRODUCT NUMBER	DESCRIPTION	LINES OF INFO	BANDS & FORMS PER BOX
2236-00-PDJ	Securlink Blood System	2	250 bands & forms
2136-00-PDJ	Securlink Blood Band (only) Red Information area: 4" L x 5/16" H		250 bands
2200-00-PDJ	Securlink Form (only)		250 forms

XTENDER BAND™

- Xtender Band™ snaps onto most PDC wristbands to offer 4-6 more inches of length for optimum adjustability and comfort
- Durable, water-resistant material
- Non-transferable, secure plastic snap provides maximum security



PRODUCT NUMBER	SIZE OF BAND	INFORMATION AREA	BANDS PER BOX
5062-11-PDJ	10" L	N/A	250

Available in White ONLY




REDUCE ERRORS AND IMPROVE PATIENT SAFETY WITH ISBT COMPLIANT LABELS

PDC can help you meet the ISBT Blood Bank standards with our complete line of blank and pre-printed labels.



[See our full selection of products on pages 104-105 of PDC's Wristband and Label Systems Catalog.](#)



GET HIGH IMPACT COMMUNICATION WITH HIGHLY VISIBLE MEDVISION® LABELS

Choose from a variety of preprinted fluorescent laser labels to communicate important alerts. Customize your own label to meet your specific needs.



[See our full selection of products on pages 121-151 of PDC's Wristband and Label Systems Catalog.](#)



ABO & Rh CONFIRMED

Date/Tech. _____

MV04FP4641

IRRADIATED

FRESH FROZEN PLASMA

DATE THAWED _____

TIME THAWED _____

EXPIRATION DATE _____

EXPIRATION TIME _____

MV04SY4662

CONTACT US TODAY!

800-435-4242

www.pdcorp.com/bloodbands

- REQUEST SAMPLES
- GET A QUOTE
- PLACE AN ORDER



VALENCIA, CALIFORNIA
 25167 Anza Drive | Valencia, Ca 91355
 Tel 800.435.4242 or 661.257.0233 | Fax 800.321.4409
 stjohnc.com | E-mail cs@stjohninc.com

SAN FERNANDO, CALIFORNIA
 13880 Del Sur Street | San Fernando, Ca 91340
 Tel 800.772.1122 or 818.897.1111 | Fax 818.686.9314
 pdcorp.com | E-mail info@pdcorp.com

EUROPE, MIDDLE EAST & AFRICA
 Rue de l'Industrie 17 | 1400 Nivelles, Belgium
 Tel +32 67 89 56 56 | Fax +32 67 21 77 70
 pdc-europe.com | E-mail info@pdc-europe.com

LATIN AMERICA & ASIA
 Tel +1.818.897.1111 | Fax +1 818.686.9317
 E-mail intl@pdcorp.com

QUALITY PROCESSES. QUALITY CARE.

As an ISO-9001 certified company, PDC follows a systematic world-standard approach to ensure superior product design, manufacturing, and customer support services. Every product we produce is born of thousands of hours of rigorous research and testing to provide the most durable, reliable, and comfortable products available. We are constantly working and innovating to ensure optimal performance and enhancements for every product we manufacture.