

# Securline® BLOOD BAND *plus*



*With over 50 years of experience, Precision Dynamics Corporation (PDC) is the leader in patient wristband identification, delivering innovative and reliable solutions to healthcare organizations worldwide.*

Providing maximum security, accuracy, and comfort, **Securline® Blood Band Plus** is the preferred choice of healthcare facilities and laboratories.

- Helps streamline processes, reduce labor, save time, and prevent processing errors.
- Meets the American Hospital Association guidelines for a tamperproof and non-transferable wristband that minimizes the risk of lost or transferred data.
- Metal-free wristband follows the recommended practice of the Association of periOperative Registered Nurses (AORN), encouraging that patients should not wear metal in the event of electrosurgery to avoid the possibility of a burn.



[www.pdcorp.com/healthcare](http://www.pdcorp.com/healthcare)



## Securline® Blood Band Plus

### 1 SecurSnap® Closure

PDC's unique SecurSnap closure provides maximum security and easily snaps together (even when personnel wear gloves), providing greater ease of use over other blood band styles.

### 2 Easy-to-Read Data Area

Information on the Blood Band is easier to read than other blood band styles, reducing the chance of a medical error. Patient data can be handwritten on the clear white label area, which is larger than competitive wristbands. Once the label is removed, a "carbonless" copy of each patient's information remains on the band.

### 3 Efficient Easy-Peel Labels

23 easy-peel labels allow medical personnel to quickly track and identify specimens and crossmatch them to the correct patient. Securline Blood Band has more ID labels than any other style of blood band. Labels conveniently remain on the wristband for use with paperwork and medical charts.

### 4 Adhesive Peel-off Tail

Tail has peel-off self-adhesive feature that applies directly to specimen tubes.

### Convenient Packaging

Securline Blood Bands are packaged in cases of 150 (15 easy-to-use dispenser boxes of 10 wristbands). Each box includes an extender band in the event the wristband needs to be applied to a patient's leg or to a larger sized person.



1 SecurSnap® Closure



2 Easy-to-Read Data Area

3 Easy-Peel Labels



4 Adhesive Peel-off Tail



Call Customer Care today for free wristband samples  
**800.772.1122**  
 6:30am to 4:30pm PST



13880 Del Sur Street  
 San Fernando, CA 91340  
 T: 818.897.1111 | F: 818.686.9314  
 info@pdcorp.com | www.pdcorp.com

### Ordering Guide

**NEW! SECURLINE BLOOD BAND PLUS** – with adhesive peel-off tail

Product Number	Name	Bands/Case	Box Weight
2406-**-PDB	Securline Blood Band Plus	150	5.70 lbs.
Substitute ** for numbered color: 16-Red, 14-Yellow, 22-Green			

**SECURLINE BLOOD BAND**

Product Number	Name	Bands/Case	Box Weight
2403-**-PDB	Securline Blood Band	150	5.70 lbs.
Substitute ** for numbered color: 16-Red, 14-Yellow, 22-Green			