

ADMISSIONS

The right label for the right application

PDC's expertise in labeling means that you don't have to waste valuable time "experimenting" with products that might not work for your application. We provide stock products, not stock answers. With over 2,000 stock products, we offer the largest selection of stock labeling systems with a 99% fulfillment rate!

Your PDC representative will find the right solution to a labeling issue in our vast inventory of stock products. Our economical products are ready for almost every application. No one else in labeling or forms boasts such a large inventory of stock labeling systems or matches our 99% fulfillment rate!



- **GENERAL PURPOSE ADMISSION LABELS**
- **HIPAA COMPLIANT LABELING SOLUTIONS**
- **FILE LABEL SOLUTIONS**
- **MEDVISION® LABELS**

ADMISSIONS PRODUCTS

Meet these Joint Commission Patient Safety Goals & Standards

Goal #1:

Improve the accuracy of patient identification.

Goal #2:

Improve the effectiveness of communication among caregivers.

Goal #8:

Accurately and completely reconcile medications across the continuum of care.



All PDC Admissions products are latex-free and some are available with custom imprinting.


GENERAL PURPOSE ADMISSION LABELS

ADMISSION LABELS





General Purpose Admission Labeling

Patient Identification labels for addressograph and other patient name plate imprinters are stocked and available in two blank white paper stocks - a smudgeproof litho for clean flat surfaces and a flexo stock for curved, irregular surfaces. Both have a removable adhesive and perforations between each label. Color-tinted ID labels or those with a permanent adhesive are available on a custom basis.






Product Number	Label Size	Core Size	Material/Adhesive	Qty.
ADMIT-1	15/16" x 3"	Fanfold	Direct Thermal/Perm	7M/BX
ADMIT-1A	15/16" x 3"	3"	Direct Thermal/Perm	5.5M/RL-1RL/BX

* Partial varnish strip helps to protect sensitive patient information.



Product Number	Label Size	Core Size	Max OD	Material/Adhesive	Qty.
ADLR-1123	1-1/2" x 3"	1"	3-7/8"	Smudgeproof Litho/Perm	1M/RL
ADFR-1123	1-1/2" x 3"		3-7/8"	Smudgeproof Flexo/Perm	
ADLR-1143	1-1/4" x 3"		3-3/4"	Smudgeproof Litho/Perm	
ADFR-1143	1-1/4" x 3"		3-3/4"	Smudgeproof Flexo/Perm	
ADLR-11432	1-1/4" x 3-1/2"		3-5/8"	Smudgeproof Litho/Perm	
ADFR-11432	1-1/4" x 3-1/2"		3-5/8"	Smudgeproof Flexo/Perm	





Product Number	Label Size	Core Size	Max OD	Material/Adhesive	Qty.
ADLR-1132	1" x 3-1/2"	1"	3-3/8"	Smudgeproof Litho/Perm	1M/RL
ADFR-1132	1" x 3-1/2"			Smudgeproof Flexo/Perm	
ADLR-13	1" x 3"			Smudgeproof Litho/Perm	
ADFR-13	1" x 3"			Smudgeproof Flexo/Perm	
ADLR-11232	1-1/2" x 3-1/2"			Smudgeproof Litho/Perm	
ADFR-11232	1-1/2" x 3-1/2"			Smudgeproof Flexo/Perm	

GENERAL PURPOSE ADMISSION LABELS







ADMISSION LABELS



General Purpose Admission Labeling



Product Number	Label Size	Core Size	Max OD	Material/Adhesive	Qty.
TTR-13	1" x 3"	3"	8"	Thermal Transfer/Glove	5.5M/RL-2RLS/BX
TTR-13Z	1" x 3"		5-9/16"	Thermal Transfer/Glove	2.5M/RL-4RLS/BX
TTR1-13	1" x 3"	1"	5"	Thermal Transfer/XPermanent	2.5M/RL-8RLS/BX
TTRS-13	1" x 3"	3"	8-3/8"	Thermal Transfer Film/Glove	5,330/RL-2RLS/BX
TTR-24	2" x 4"		7-1/2"	Thermal Transfer/Glove	1M/RL-6RLS/BX



Product Number	Label Size	Core Size	Max OD	Material/Adhesive	Qty.
PSLP-12	1" x 2"	1"	3-1/2"	Smudgeproof Litho/Perm	1M/RL
PSLP-1212	1" x 2-1/2"		3-3/8"		
PSLP-1214	1" x 2-1/4"		3-1/2"		
PSLP-11212	1-1/2" x 2-1/2"		4"		
PSLP-13424	1-3/4" x 2-3/4"		4-1/4"		
PSLP-23	2" x 3"		4-1/2"		

COMPANION ADMISSION LABELING SOLUTIONS

ADMISSION LABELS



Kits: Label and Ribbon Combinations



These labels are ideal for applications that require durability or the highest level of resistance to abrasion, heat, steam and various solvents. Thermal labels are often used for Laboratory and Central Service/Sterile applications.


Kits: Label and Ribbon Combinations					
Product Number	Label/Ribbon Size	Core Size	Max OD	Material/Adhesive	Qty.
TTR-1583KIT	1-5/8" x 3"	3"	7-1/2"	Thermal Transfer/Glove	3M/RL-3RLS/BX
	3.15" x 1,476'	-	-	Resin Enhanced Wax	1/RL/BX

See the Thermal Transfer Product section of this catalog for more information

HIPAA Compliant Labeling Solutions



These unique block out labels have a "black" adhesive which blocks data from being read from the front and the back. The permanent black adhesive also destroys patient data upon removal.



Product Number	Label Size	Core Size	Material/Adhesive	Qty.
MV02SW8216	1-7/16" x 3/8"	1"	Smudgeproof Litho/Perm	666/RL
MV04SW8215	2-1/4" x 1"			420/RL
MV06SW8217	2-15/16" x 1"			333/RL
MV08SW8218	2-15/16" x 2"			333/RL
MV09SW8219	2-5/8" x 4"			375/RL

Visitor Pass Labels

For the safety of your facility, securely record visitor traffic by issuing pre-printed visitor pass labels.



Product Number	Label Size	Core Size	Material/Adhesive	Qty.
TM-VP-RED	1-3/4" x 2-3/4"	1"	Name badge/Rem	500/RL-2RLS/BX



Basic File Labels - Product Mix

Please keep in mind that this is just a sample product mix. The full product line and product details can be found in the Filing System Labels brochure available through our Customer Care department or at www.timemed.com.

Alpha Labels

Prefix Definitions

AL = Alpha Labels
WR = Write-In Boxes

NL = Numeric Labels
YR = Year Labels

SL = Solid Labels

1. Alphabetical Numbering System -

Ideal for less than 15,000 patient charts. Files by the individual's last name with no additional cross-referencing.

Pros - Easiest to use since all records are filed by patient's last name.

Cons - Exact placement of records is required so that alpha sequence is maintained.



Numeric Labels

2. Straight Numeric Numbering System -

Ideal for fewer than 10,000 patient charts. Files in sequential order starting from the lowest to highest number; alphanumeric cross-reference is required.

Pros - Arranged from lowest to highest number - provides easy way to retrieve records that are assigned in consecutive order.

Cons - Requires computerized cross-reference to patient name in order to locate record.

3. Terminal Digit Numbering System -


Ideal for approximately 100,000-200,000 patient charts. Files in groups of like numbers segmenting the original number into smaller groups of 100; alphabetical cross-reference is required.

Pros - Activity is equally allocated throughout the system. This eliminates concentrated activity and provides equal expansion and reduction (purging) within the system.

Cons - Requires computerized cross-reference to patient name in order to locate record.




Date Labels



Write-In Boxes

Use these labels to cover digits so you can reuse folders with a new number. One label for each side of folder tab.

Available in
1/4" and 1/2" sizes



TimeMed MedVision® Possibilities

ALLERGY

MV04FR001

Protected Health Information
DO NOT RELEASE

MV06FR213

MEDVISION LABELS FOR ADMISSIONS

ADMISSION LABELS



The MedVision® Label Program provides a quick and convenient way to order a wide variety of preprinted or custom labels for Admissions. The labels shown on this page represent a small selection of the many MedVision® labels available. View hundreds of stock labels in the MedVision® section (page 121) or view the complete selection at www.timemed.com.

NAME ALERT
MV01FR1207

NO CAFFEINE
MV01FR0531

Food Allergies

MV03FR8763

ATTENTION

HEARING IMPAIRED

SPEECH IMPAIRED

VISION IMPAIRED

FOREIGN LANGUAGE

MV03FP5467

Product Number	Description	Label Size	Color	Material/Adhesive	Qty.
MV01FR1207	Name Alert	1-1/4" x 5/16"	Fluorescent Red	Fluorescent Red/Permanent	760/RL
MV01FR0531	No Caffeine	1-1/4" x 5/16"	Fluorescent Red	Fluorescent Red/Permanent	760/RL
MV03FR8763	Food Allergies	1-7/16" x 1"	Fluorescent Red	Fluorescent Red/Permanent	666/RL
MV03FP5467	Attention	1-7/16" x 1"	Fluorescent Pink	Litho/Removable	666/RL

CHART THINNED

On: _____

By: _____

MV04FG1212

**Fluid warming cabinet has been used.
Do NOT return to warming cabinet
after removal. Remove from warmer by:**

Date: _____ Initials: _____


MV04FO5924

POOR ORIGINAL DOCUMENT

MV04SW6885

Product Number	Description	Label Size	Color	Material/Adhesive	Qty.
MV04FG1212	Chart Thinned	2-1/4" x 1"	Fluorescent Green	Litho/Removable	420/RL
MV04FO5924	Fluid Warming		Fluorescent Orange	Litho/Removable	
MV04SW6885	Poor Original		White	Litho/Removable	

DOCTOR PLEASE NOTE



DATE _____

MV05FP7051

**Maintain as a part of permanent
medical record.**

DATE REPORTED: _____

TIME REPORTED: _____

PERFORMED BY: _____

REPORTED TO: _____
(Name and Credential)

MV07SW8842

Product Number	Description	Label Size	Color	Material/Adhesive	Qty.
MV05FP7051	Doctor Please Note	2-7/16" x 2-1/2"	Fluorescent Pink	Litho/Removable	400/RL
MV07SW8842	Maintain as part of	2-15/16" x 1-1/2"	White		333/RL



TO: _____
FROM: _____
PLEASE NOTE

_____ **DATE:** _____

MESSAGE NOTED: _____

Medicare

Expected date of discharge _____

1st date IM delivered _____

2nd date IM delivered _____

Appeal date (if applicable) _____

MV09SY0766

Product Number	Description	Label Size	Color	Material/Adhesive	Qty.
MV09SW0218	To: From: Please Note	2-5/8" x 4"	White	Litho/Removable	375/RL
MV09SY0766	Medicare		Yellow	Litho/Permanent	

Our 1052P-B dispenser provides the perfect tool to keep your MedVision® labels organized

PDC dispensers are designed to provide convenient dispensing of any tape or label. Manufactured from sturdy plastic, they will provide a lifetime of use. Set them on your shelf, mount them on a wall, hang them from your crash cart or hold them in your hand!



Product Number	Description
1052P-B	Can hold multiple rolls of MedVision® labels up to a total of 8-1/2". Includes inserts, hanging bars and bottom suction cups for secure placement. Size: (H) 4-5/8" x (W) 9" x (D) 6"